

Recovery:

Known amounts of AMH/MIS were added to three serum samples containing different levels of endogenous AMH/MIS. The concentration of AMH/MIS was determined before and after the addition of exogenous AMH/MIS and the percent recovery was calculated.

Sample	Endogenous Conc. (ng/mL)	Expected Conc. (ng/mL)	Observed Conc. (ng/mL)	% Recovery
1	1.8408	2.329	2.257	97%
		2.817	2.598	92%
		3.305	3.039	92%
2	0.3218	0.886	0.821	93%
		1.450	1.250	86%
		2.014	1.749	87%
3	1.2683	1.785	1.676	94%
		2.301	2.027	88%
		2.818	2.493	88%

ANALYTICAL SPECIFICITY

This monoclonal antibody pair used in the assay detects human AMH/MIS. Specificity to other species has not been determined.

Cross-Reactant	Concentration	% Cross-reactivity
Inhibin A	100 ng/mL	ND
Inhibin B	100 ng/mL	ND
Activin A	50 ng/mL	ND
Activin B	50 ng/mL	ND
Activin AB	50 ng/mL	ND
rAMH	130 ng/mL	ND
Mature AMH	120 ng/mL	2.35
hAMH(Pro)	300 ng/mL	0.24
ProMature hAMH	110 ng/mL	0.09

Interference:


When potential interferents (hemoglobin, triglycerides and bilirubin) were added at least at two times their physiological concentration to control sample, AMH/MIS concentration were within $\pm 10\%$ of the control as represented in the following table. This study was based on NCCLS EP7-P to serum matrix added.

Interferents	Analyte Conc. (mg/mL)	Unspiked Sample Value (ng/mL)	Spiked Sample Value (ng/mL)	% Difference
Hemoglobin	1.35	5.9979	5.8795	-1.974
		9.0887	9.081	-0.085
Triglycerides	5.0	5.9979	6.0296	0.529
		9.0887	8.8229	-2.925
Bilirubin	0.6	0.7818	0.8009	2.443
		9.3466	9.0302	-3.385

Sample Type:


Forty matched serum and Lithium heparin plasma specimens in the range of 0.14-15.14 ng/mL were compared in AnshLite™ US AMH/MIS CLIA assay. Passing Bablok analysis of the results yielded the following Regression:

$$\text{Plasma} = 1.03 (\text{serum}) - 0.12, (r=0.995; P<0.0001)$$

**Method Comparison:**

The UltraSensitive AMH/MIS CLIA has been compared to Commercial AMH assay (Method A) using 90 serum samples. Passing Bablok analysis of the results yielded the following Regression:

$$\text{US AMH/MIS CLIA (AL405)} = 1.10 (\text{Method A}) + 0.10 \\ (r=0.915; P<0.0001)$$

**REFERENCES**

1. HHS Publication, 5th ed., 2007. Biosafety in Microbiological and Biomedical Laboratories. Available <http://www.cdc.gov/biosafety/publications/bmbL5/BMBL5>
2. DHHS (NIOSH) Publication No. 78-127, August 1976. Current Intelligence Bulletin 13 - Explosive Azide Hazard. Available <http://www.cdc.gov/niosh>.
3. Approved Guideline – Procedures for the Handling and Processing of Blood Specimens, H18-A3. 2004. Clinical and Laboratory Standards Institute.
4. Kricka L. Interferences in immunoassays – still a threat. Clin Chem 2000; 46: 1037-1038.

FOR RESEARCH USE ONLY

Not for use in diagnostic procedures.

The Ansh Labs logo is a trademark of Ansh Labs.



Manufactured by:

Ansh Labs

445 Medical Center Blvd.

Webster, TX 77598-4217, U.S.A.