

7. During the last **20-30 minutes** of incubation, prepare the Inhibin B Antibody-Biotin Conjugate Solution by diluting the Inhibin B Biotin Conjugate Concentrate in Inhibin B Conjugate Diluent as described under the Preparation of the Reagents section of this package insert.
8. Aspirate and wash each strip **5 times** with Washing Solution (**350 µL/per well**) using an automatic microplate washer.
9. Add **100 µL** of the Antibody-Biotin Conjugate Solution to each well using a repeater pipette.
10. Incubate the plate, shaking at a fast speed (**600-800 rpm**) on an orbital microplate shaker, for **1 hour** at room temperature.
11. Aspirate and wash each strip **5 times** with the Wash Solution (**350 µL/per well**) using an automatic microplate washer.
12. Add **100 µL** of the Streptavidin-Enzyme Conjugate-RTU to each well using a repeater pipette.
13. Incubate the plate, shaking at a fast speed (**600-800 rpm**) on an orbital microplate shaker, for **30 minutes** at room temperature.
14. Aspirate and wash each strip **5 times** with the Wash Solution (**350 µL/per well**) using an automatic microplate washer.
15. Add **100 µL** of the Substrate Solution to each well using a repeater pipette.
16. Incubate the wells, shaking at **600-800 rpm** on an orbital microplate shaker for 4±1 minute at room temperature. Avoid exposure to direct sunlight.
17. Read the light output of the solution in the wells within 10 minutes, using a microplate luminometer.
NOTE: While reading the RLU'S of the well, it is necessary to program the 0 calibrator as "BLANK."

ALTERNATIVE ASSAY PROCEDURE

1. Reconstitute Inhibin B Calibrator A-F and Inhibin B Controls I & II each with **1 mL** deionized water. Solubilize for **10 minutes**, Mix well.
2. For each unknown serum sample, calibrators and controls label one 12 X 75 culture tubes.
3. Pipette **75 µL** of the Calibrator, Controls and samples to the pre-labeled tube.
4. Add **75 µL** of the Inhibin B Assay Buffer A to each the pre-labeled tube using a repeater pipette.
5. Add **75 µL** of the Inhibin B Assay Buffer B to each pre-labeled tube using a repeater pipette and vortex well.
6. Place the tubes in a tight fitting tube rack and incubate the tubes, shaking at a slow speed (**100-200 rpm**) at room temperature for **30 minutes**.
7. The pre-mixed samples are now ready for analysis.
8. Label the microtitration strips to be used.
9. Pipette **150 µL** of the pre-mixed samples from **step 7** to the appropriate wells.
10. Incubate the plate, shaking at a fast speed (**600-800 rpm**) on an orbital microplate shaker, for **2 hour** at room temperature.
11. During the last **20-30 minutes** of incubation, prepare the Inhibin B Antibody-Biotin Conjugate Solution by diluting the Inhibin B Biotin Conjugate Concentrate in Inhibin B Conjugate Diluent as described under the Preparation of the Reagents section of this package insert.
12. Aspirate and wash each strip **5 times** with Washing Solution (**350 µL/per well**) using an automatic microplate washer.
13. Add **100 µL** of the Antibody-Biotin Conjugate Solution to each well using a repeater pipette.
14. Incubate the plate, shaking at a fast speed (**600-800 rpm**) on an orbital microplate shaker, for **1 hour** at room temperature.
15. Aspirate and wash each strip **5 times** with the Wash Solution (**350 µL/per well**) using an automatic microplate washer.
16. Add **100 µL** of the Streptavidin-Enzyme Conjugate-RTU to each well using a repeater pipette.

17. Incubate the plate, shaking at a fast speed (**600-800 rpm**) on an orbital microplate shaker, for **30 minutes** at room temperature.
18. Aspirate and wash each strip **5 times** with the Wash Solution (**350 µL/per well**) using an automatic microplate washer.
19. Add **100 µL** of the Substrate Solution to each well using a repeater pipette.
20. Incubate the wells, shaking at **600-800 rpm** on an orbital microplate shaker for 4±1 minute at room temperature. Avoid exposure to direct sunlight.
21. Read the light output of the solution in the wells within 10 minutes, using a microplate luminometer.
NOTE: While reading the RLU'S of the well, it is necessary to program the 0 calibrator as "BLANK."

RESULTS

NOTE: The results in this package insert were calculated by plotting the data on a log vs. log scale using a cubic regression curve-fit. Other data reduction methods may give slightly different results.

1. Calculate the mean RLU for each calibrator, Control, or Unknown.
2. Plot the log of the mean RLU readings for each of the Calibrators along the y-axis versus log of the Inhibin B concentrations in pg/mL along the x-axis, using a cubic regression curve-fit.
3. Determine the Inhibin B concentrations of the Controls and unknowns from the calibration curve by matching their mean RLU readings with the corresponding Inhibin B concentrations.
4. Any sample reading higher than the highest Calibrator should be appropriately diluted with the 0 pg/mL (CAL A) and re-assayed.
5. Any sample reading lower than the analytical sensitivity should be reported as such.
6. Multiply the value by a dilution factor, if required.

LIMITATIONS

The reagents supplied in this kit are optimized to measure Inhibin B levels in human serum. If there is evidence of microbial contamination or excessive turbidity in a reagent, discard the vial. For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the samples⁴.

QUALITY CONTROL

- Each laboratory should establish mean values and acceptable ranges to assure proper performance.
- Inhibin B CLIA controls or other commercial controls should fall within established confidence limits.
- The confidence limits for Inhibin B CLIA controls are printed on the **Calibration card**.
- A full calibration curve, low and high level controls, should be included in each assay.
- The Substrate solutions should be colorless. Development of any color may indicate reagent contamination or instability.

REPRESENTATIVE CALIBRATION CURVE DATA

Well Number	Well Contents	Mean RLU	Conc (pg/mL)
Calibrators			
A1, A2	A	0.318 (Blank)	0
B1, B2	B	0.437	14
C1, C2	C	1.993	42
D1, D2	D	9.638	125
E1, E2	E	45.841	415
F1, F2	F	158.287	1344

CAUTION: The above data must not be employed in lieu of data obtained by the user in the laboratory.

