

DESCRIPTION

Specificity	GDF-15, H202 mutation
-------------	-----------------------

Source	Murine
--------	--------

Product Form / Appearance	Clear; Colorless solution
---------------------------	---------------------------

Product Type	Monoclonal Antibody
--------------	---------------------

Isotype	IgG2a
---------	-------

Package Size	100 µg
--------------	--------

SPECIFICATIONS

Species Reactivity	Human, mouse
--------------------	--------------

Immunogen	Recombinant GDF-15
-----------	--------------------

Purification	Purified from hybridoma culture supernatant using Protein G.
--------------	--

Concentration	Batch Specific – refer to Certificate of Analysis report
---------------	--

Activity	Determined by direct ELISA
----------	----------------------------

Composition	Reconstitute to the appropriate working dilution for your intended application using 10mM Phosphate Buffered Saline at pH 7.4 without preservatives.
-------------	--

PREPARATION AND STORAGE

Reconstitution	Not Applicable. Provided in liquid format.
----------------	--

Shipping Conditions	This product is shipped on dry ice. Upon receipt, store it immediately at the temperature recommended below.
---------------------	--

Storage	Store at -20°C to -80°C as supplied. Avoid repeated freeze-thaw cycles.
---------	---

APPLICATIONS

This product has been reported to work in the following applications. This information is derived from testing within our laboratory, peer-reviewed publications or personal communications from the originators. Optimal dilutions should be determined by each laboratory for each application.

	Status	Recommended Concentration
ELISA	Recommended	10 µg/mL
SDS-Page	Recommended	1-10 µg