

Pricing: **\$800 - 5µg***Inquire about other package sizes and bulk pricing.***DESCRIPTION**

Specificity	PAPP - A
Source	Human Placental Fluid
Product Form / Appearance	Clear liquid solution
Product Type	Antigen
Application	ELISA

**SPECIFICATIONS**

Concentration	Batch Specific - Refer to C of A
Activity	Determined by ELISA
Composition	10mM TRIS buffer pH7.4, 0.15M NaCl and 0.09% Azide
Package size	5 µg

**STORAGE, HANDLING, HAZARDS**

Reconstitution	Reconstitution not applicable, product supplied as liquid.
Shipping Conditions	This product is shipped on dry ice. Upon receipt, store it immediately at the temperature recommended below.
Storage	Store at -20°C to -80°C as supplied.
Chemical hazard	This product does not contain any hazardous chemicals.
Biological hazard	This product is not known to pose any infectious risk. However, no available test method can offer complete assurance of eliminating potential biohazardous risk and hence, this product should be handled with care as a "hazardous component."