

Pricing: **\$800 - 5µg***Inquire about other package sizes and bulk pricing.***DESCRIPTION**

|                           |                      |
|---------------------------|----------------------|
| Specificity               | PAPP - A2            |
| Source                    | CHO-S Cells          |
| Product Form / Appearance | Pink liquid solution |
| Product Type              | Antigen              |
| Application               | ELISA                |

**SPECIFICATIONS**

|               |                                  |
|---------------|----------------------------------|
| Concentration | Batch Specific - Refer to C of A |
| Activity      | Determined by ELISA              |
| Composition   | Supernatant, No Preservatives    |
| Package size  | 5 µg                             |

**STORAGE, HANDLING, HAZARDS**

|                     |  |
|---------------------|--|
| Reconstitution      | Reconstitution not applicable, product supplied as liquid.   |
| Shipping Conditions | This product is shipped on dry ice. Upon receipt, store it immediately at the temperature recommended below.   |
| Storage             | Store at -20°C to -80°C as supplied.   |
| Chemical hazard     | This product does not contain any hazardous chemicals.   |
| Biological hazard   | This product is not known to pose any infectious risk. However, no available test method can offer complete assurance of eliminating potential biohazardous risk and hence, this product should be handled with care as a "hazardous component." |